# PUNKTAL

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**Digital**

**PK-43BAL**

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## INSTRUCTIONS FOR USE / ENGLISH

### KEY FEATURES

1. High precision strain gauge sensors system .
2. Tempered-glass platform.
3. Capacity: 6kg-180kg (Unit:Kg Lb St:lb)
4. Graduation: d= 100g
5. Auto on
6. Auto zero & Auto off
7. Low battery and over load indication
9. Free switch three different units: kg/lb/st:lb.
10. Powered by 1\*3V lithium cell

### OPERATION

1. Take the insulating strip out

2. Place the scale on hard and flat surface. Clean floor is best.

3. Step on the scale gently, position your feet evenly near the center and keep your body still and erect until the scale stabilizes your weight. The scale will be auto off within seconds if no further using.

### CAUTION

- To ensure accuracy, we suggest you take measurement at the same time every night (19:00-21:00).
- Please change new batteries when “LO” appears on the LCD which indicates that the battery is out of power.

- “Err” showed on LCD indicates overload and “Err” indicates that the measuring makes an error.

### ATTENTION

- Please take off shoes before taking measurement to avoid scraping or scratching the glass platform.
- The Glass platform may be pretty slippery when wet, thus please make sure both the glass your feet are dry before using.
- Please use soft tissue with alcohol or glass cleaner to clean the surface if it is dirty. No soap or other chemicals. Keep it away from water, heat and extreme coldness.
- The scale is high precision measuring device, Never jump or stomp on the scale or disassemble it and please handle it carefully to avoid dropping when moving it.
- The scale is for family use only and not suitable for professional use.
- Please check the battery if the scale malfunctions. Chang a new piece if necessary

IMPORTADOR EN URUGUAY:  
**PUNKTAL S.A.**  
RONDEAU 1999 – MONTEVIDEO  
TEL: 2929 1717  
[www.punktal.com.uy](http://www.punktal.com.uy)

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